

NEW for the era of liquid automation...

NOW AVAILABLE

Probact Liquid Amies Transport Swabs
for integration with high volume automated laboratories



- > Easy collection of aerobes, anaerobes and fastidious organisms
- > Flocked swab for reliable absorbency and release of sample
- > Leak proof, swab capture cap for increased security
- > Compatible with most automated laboratory equipment
- > Suitable for both culture & molecular assays
- > CE marked and M40 compliant

SAMPLES AVAILABLE FOR EVALUATION



Technical Service
Consultants Ltd

www.tscswabs.co.uk

Probact Liquid Amies Transport Swab



Probact Liquid Amies Transport Kits are for the collection and transport of bacterial samples for integration with high volume automated laboratories.

CE marked and CLSI M40-A compliant, the all in one sampling kit contains 1ml of liquid Amies medium, ideal for maintaining the viability of aerobic, anaerobic, and fastidious bacteria for up to 48 hours and a superior flocked swab for the reliable collection and transport of bacteria to the laboratory.

Simple to use

Swab kit includes a vial with 1ml of liquid Amies medium and flocked swab in a tamper evident peel pouch.

Sample the surface

Using the 15cm long polystyrene shaft with soft flocked swab tip, sample absorption and release is optimised by the varied length nylon fibres and is ideal for general sampling including nasopharyngeal, endocervical and wound sampling.

Immerse the swab

Place the flocked swab into the liquid medium in the vial.

Snap the swab

Bend the swab in the vial and it will conveniently snap at its breakpoint.

Secure the swab

Place the swab capture cap on the vial and tighten, securing the swab to the cap.

Send to the laboratory

Samples should be sent for analysis as soon as possible. Viability of aerobic, anaerobic, and fastidious bacteria will be maintained for up to 48 hours (except for *Neisseria gonorrhoeae*, which should be processed within 24 hours).

Code	Product	Description	Pack Size
TS/6-17	Probact Liquid Amies Transport Swab	1ml Liquid Amies media in vial and flocked swab with flexible polystyrene shaft	50
HHD298 NHS NPC Code	Probact Liquid Amies Transport Swab	1ml Liquid Amies media in vial and flocked swab with flexible polystyrene shaft	50

ISO 9001:2008

EN ISO 13485:2012

CLSI M40-A Compliant

Medical Device Directive 93/42/EEC



**Technical Service
Consultants Ltd**

Technical Service Consultants Ltd

Microbiology House, Fir Street, Heywood, Lancashire OL10 1NW, UK

T : +44 (0)1706 620600 F : +44 (0)1706 620445

E : sales@tscswabs.co.uk W : www.tscswabs.co.uk