

Performance Characteristics & References

SwabSURE ListeriaP

Performance Characteristics

Refer to the following table for examples of SwabSURE-ListeriaP Kit performance data. Detailed specificity, sensitivity and swab performance data are available on request.

Organism	Test reaction	Colour reaction	Concentration, CFU/Swab
<i>Listeria monocytogenes</i>	Positive	Turquoise-blue	≥ 1
<i>Listeria ivanovii</i>	Positive	Turquoise-blue	≥ 1
<i>Listeria seeligeri</i> ⁴	Negative/Positive	Yellow/Turquoise-blue	10-1000
<i>Listeria innocua</i>	Negative	Yellow	10-1000
<i>Listeria welshimeri</i>	Negative	Yellow	10-1000
<i>Bacillus cereus</i>	Negative	Yellow	10-1000
Lactobacillus spp	Negative	Yellow	10-1000
<i>Staphylococcus aureus</i>	Negative	Yellow	10-1000
Micrococcus spp	Negative	Yellow	10-1000
Enterobacter spp	Negative	Yellow	10-1000
Enterococcus spp	Negative	Yellow	10-1000
Candida spp	Negative	Yellow	10-1000

References

1. Manafi M. 2000. New developments of chromogenic and Fluorogenic culture media. Int. J. Food Microbiol. 60: 205-218.
2. Notermans S.H.W., Dufrenne J., Leimeister-Wächter M., Domann E., Chakraborty T. 1991. Phosphatidylinositol-specific phospholipase C activity as a marker to distinguish between pathogenic and nonpathogenic *Listeria* species. Applied and Environmental Microbiology. 57 (9): 2666-2670.
3. Restaino L., Frampton E.W., Irbe R.M., Schabert G., Spitz H. 1999. Isolation and detection of *Listeria monocytogenes* using fluorogenic and chromogenic substrate for phosphatidylinositol-specific phospholipase C. J. Food Prot. 62, 244-51.
4. Rocourt, J., H. Hof, A. Schrettenbrunner, R. Malinverni, and J. Bille. 1986. Acute purulent *Listeria seeligeri* meningitis in an immunocompetent adult. Schweiz Med. Wochenschr. 116:248-251.

Have you considered? Sterile Sampling Template - TS/15-T



Technical Service
Consultants Ltd

www.tscswabs.co.uk
Tel: +44 (0)1706 620600

PC+R_v1 – 24/06/2014