



Product Specification - **TS/12-A**, 1 µl Blue Culture Loops (soft), 20 in the sealed minigrip bag

Product Code	TS/12-A
Description	1 µl Blue Culture Loops (soft), 20 in the sealed minigrip bag
Intended Use	Sample throughout. General laboratory use only
Components	1 µl Plastic Blue Loops (soft) Minigrip bag
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	1000 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 09.04.2008

