



Product Specification - **TS/15-B: NaCl**, Product Name - 50cm² sponge with mini grip bag, individually packed & dosed with 10 ml of unbuffered saline (sodium chloride) solution

Product Code	TS/15-B: NaCl
Description	Product Name - 50cm ² sponge with mini grip bag, individually packed & dosed with 10 ml of unbuffered saline (sodium chloride) solution
Intended Use	Sampling of meat carcasses for microbiological analysis
Components	Cellulose blue sponge 50cm Mini grip bag Medical grade plastic bag
Contents	10 ml of 0.9% NaCl Sodium Chloride Solution:De-ionised Water – 1 Litre NaCl – 9 g
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	12 months from date of manufacture Storage Conditions - 4-25 0C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100
Other	Quality control testing details on final product: i) Validation of irradiation sterility.ii) Product Specification iii) Organism recovery

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 29.04.2008

