



Product Specification - **TS/15-B: PSD**, Sterile 50cm<sup>2</sup> blue sponge with mini grip bag, individually packed & dosed with 10 ml Peptone Saline Diluent (PSD)

Product Code	TS/15-B: PSD
Description	Sterile 50cm <sup>2</sup> blue sponge with mini grip bag, individually packed & dosed with 10 ml Peptone Saline Diluent (PSD)
Intended Use	Sampling of meat carcasses for microbiological analysis
Components	Cellulose blue sponge 50+2cm <sup>2</sup> Mini grip bag Medical packaging grade sponge bag Label and Film Peptone Saline Diluent (MRD)
Contents	De-ionised Water – 1 Litre Bacteriological Peptone – 1g NaCl – 8.5 g
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	150
Other	Quality control testing details on final product: i) Validation of irradiation sterility.ii) Organism recovery (S.typhimurium)

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 16.10.2007

