



Product Specification - **TS/15-BPW**, Blue 50cm² sponge with mini grip bag, individually packed & dosed with Buffered Peptone Water

Product Code	TS/15-BPW
Description	Blue 50cm ² sponge with mini grip bag, individually packed & dosed with Buffered Peptone Water
Intended Use	Swab for pathogen monitoring in poultry houses
Components	Small blue sponge 50cm ² Mini grip bag Medical packaging grade sponge bag Peptone Water
Contents	De-ionised Water Bacteriological Peptone Water Powder Phosphate Buffers
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	12 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	150
Other	Quality control testing details on final product: i) Validation of irradiation sterility.ii) Organism Recovery (S.typhimurium)

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 09.04.2008

