



Product Specification - **TS/5-1**, 5ml Amies transport medium without charcoal in tube (blue cap) with blue cap cotton wool tipped(wood shaft) transport swab in a peel pouch

Product Code	TS/5-1
Description	5ml Amies transport medium without charcoal in tube (blue cap) with blue cap cotton wool tipped(wood shaft) transport swab in a peel pouch
Intended Use	This product is intended for bacteriological samples only DO NOT USE THIS PRODUCT WHEN CHLAMYDIAL INFECTION IS EXPECTED
Components	Blue cap Polypropylene Tube Blue cap cotton/wood swab Non-toxic glue Label & film
Contents	De-ionised water Sodium chloride Disodium hydrogen phosphate Sodium thioglycollate Calcium chloride Magnesium chloride Potassium phosphate Agar
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100
Other	Quality control testing details on final product: i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 29.01.2008

