



Product Specification - **TS/5-13**, 5ml Amies transport medium without charcoal in tube with blue cap cotton tipped (polystyrene shaft)swab in a peel pouch

Product Code	TS/5-13
Description	5ml Amies transport medium without charcoal in tube with blue cap cotton tipped (polystyrene shaft)swab in a peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Blue cap Polypropylene Tube Blue cap cotton/polystyrene (breakable) swab Non-toxic glue Label & film
Contents	De-ionised water Amies Transport Medium pH – 7.15 ± 0.2
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100
Other	Quality control testing details on final product: i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 15.02.2008

S:/QC/Product Specifications (inc. template)ProductSpec\_TS5-13-i1

