



Product Specification - **TS/5-18**, 5ml Amies with charcoal transport medium in tube with black cap pure viscose tipped (polystyrene shaft) swab in a peel pouch Latex free

Product Code	TS/5-18
Description	5ml Amies with charcoal transport medium in tube with black cap pure viscose tipped (polystyrene shaft) swab in a peel pouch Latex free
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Black cap Polypropylene Tube Black cap viscose/polystyrene swab Non-toxic glue Label & film Transport medium
Contents	De-ionised water Amies Charcoal Transport Medium pH – 7.15 ± 0.2
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100
Other	Quality control testing details on final product: i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 29.03.2007

S:/QC/Product Specifications (inc. template)ProductSpec_TS5-18-i2

