



Product Specification - **TS/5-18A**, 4ml Amies with charcoal transport medium in tube with black cap pure cotton tipped polystyrene shaft swab in a peel pouch

Product Code	TS/5-18A
Description	4ml Amies with charcoal transport medium in tube with black cap pure cotton tipped polystyrene shaft swab in a peel pouch
Intended Use	Bacteriological samples only. Collection and transport of clinical specimens
Components	Black cap Polypropylene Tube (80 mm) Black cap cotton tipped plastic swab (shaft length 80 mm) Non-toxic glue Label & film
Contents	De-ionised water Amies Charcoal Transport Medium pH – 7.15 ± 0.2
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Organism Recovery (N.gonorrhoeae, H.influenzae)

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 30.01.2008

