

Amies Charcoal M40 Transport Swabs

**NOW
AVAILABLE**

**In stock and available now
from the NHS Supply Chain**

Complete sampling kit for the simple and safe collection, transport and preservation of micro-organisms to the laboratory



HHD301 – Packed in 50's
HHD299 – Packed in 100's

- > Up to 48 hours survival
- > Polystyrene shaft with viscose tip
- > CE marked class II medical device
- > M40 transport media validated for excellent recovery
- > Round bottom, shatter proof tubes increase integrity

TS/6-17 - NEW Liquid Amies Transport Swab available NOW - **HHD298**



**Technical Service
Consultants Ltd**

www.tscswabs.co.uk