

The following product has been updated:

TSC Neutralising Buffer Solution

Technical Service Consultants (TSC) has revised the Neutralising Buffer pH specification to 7.2 ± 0.2 . This change affects all liquid Enviroscreen products.

This adjustment in pH contributes to improved micro-organism recovery, particularly benefiting stressed cells whilst still maintaining excellent neutralising performance.

Products updated include:

TS/5-40, TS/5-42, TS/5-43, TS/5-44, TS/5-61, TS/5-61:BP, TS/15-A:74, TS/15-B, TS/15-BG, TS/15-I, TS/15-N, TS/15-U

The change in formulation will be in effect from 17th December 2013. If you have questions regarding this change please contact us on +44 (0) 1706 620600