

Cleaning Validation Swab Kits

TSC Swab Sampling Kits contain a 15cm (6 inch) flexible swab, ideal for hard to reach areas on equipment and surfaces.

Technical Service Consultants (TSC) Swab Sampling Kits (aka Swab Rinse Kits) are all in one sampling kits for routine microbiological and general hygiene monitoring of equipment and surfaces, ideal for validation of cleaning programmes. The peel pouch kit includes a 15cm (6inch) swab shaft with choice of swab tip and a transport tube containing neutralising buffer.

Part of TSC's successful Enviroscreen Range, TS/5-42 is a sterile swab sampling kit, in a convenient and cost effective peel pouch, which is tamper evident for security assurance. The user friendly kit contains a safe and reliable 15cm (6inch) flexible polystyrene swab shaft with blue cap, to confidently reach those awkward areas on equipment and surfaces. Combined with a polyester tip, this gives excellent collection and release capabilities.

The transport tube contains 10ml of Neutralising Buffer (ideal for heavier contamination), that will neutralise commonly used disinfectants and sanitisers, allowing cleanliness to be precisely verified and maintained with minimal effort and cost.

Reach your goals confidently when testing potential contact surfaces including walls, floors, drains, utensils, containers, conveyors, whether in food and beverage production, storage and distribution, pharmaceutical and cosmetics production and on clinical contact surfaces, in-fact anywhere where you need validation of cleaning programmes.

All kits comply with ISO 18593:2004 requirements, packed in 100's for convenience and can be stored at room temperature, with a shelf life of 2 years.

For further details about the Enviroscreen Range and samples please call +44(0)1706 620600 or visit our website www.tscswabs.co.uk

