

Press Release

Release date: Feb 2012

SwabSURE

ListeriaP

The first direct detection kit for pathogenic Listeria



Technical Service Consultants SwabSURE ListeriaP is an innovative, Campden BRI validated, ISO 18593 compliant, easy to use colour-change testing system for presumptive detection of pathogenic Listeria from food contact and environmental surfaces.

ListeriaP permits easy differentiation between pathogenic *Listeria monocytogenes* and *Listeria ivanovii* from other background microorganisms including *Bacillus* spp., *Enterococcus* spp., *Micrococcus* spp., *Klebsiella* spp., and other *Listeria* species like *Listeria innocua*., simplifying the testing process involved in meeting safety and regulatory requirements such as HACCP.



The pre-moistened foam swab optimises sample collection whilst the unique TSC neutralising buffer inactivates any disinfectant residues including QAC's, phenols, peroxides and most sanitisers, enhancing cell recovery ensuring accurate and confident results every time. After the simple, safe and secure swab transfer into the **ListeriaP** medium using the convenient breakpoint, the structure of foam swab maximises sample release ensuring optimal recovery of target cells.

The selective media with its patented chromogenic substrate reliably detects *Listeria monocytogenes* down to 1 CFU by changing colour from straw to turquoise within 24 to 48 hours. This visual screen reduces unnecessary and costly subculturing and further identification needs, eliminating false positives, a common problem with other traditional methods. **ListeriaP** is simple, quick and ready to use, with the kit containing everything to perform 81 tests.



To order your Swab kit for the detection of pathogenic Listeria quote order number SS-L01. For further details, validation results or to request a free sample pack of 5 tests please visit www.tscswabs.co.uk/listeria or ring 01706 620 600.



**Technical Service
Consultants Ltd**

Technical Service Consultants Ltd

Microbiology House, Fir St, Heywood, Lancashire OL10 1NW, UK

T : +44 (0)1706 620600 F : +44 (0)1706 620445

E: sales@tscswabs.co.uk W: www.tscswabs.co.uk

11/12/2012 version8