

Technical Service Consultants Proback Range Offers Full M40 Compliance

Technical Service Consultants (TSC) Proback Transport Swab Range offers a complete sampling solution with Full M40 compliance, for the safe collection, transport and preservation of microorganisms to the laboratory.

All Proback Transport Swab Kits are Fully M40-A compliant, CE Marked and comply with the European Medical Device Directive 93/42/EEC.

- 2 year shelf life
- Tamper evident packaging
- Synthetic or natural swab materials
- High performance and cost saving options
- Excellent choice of swab heads and shaft materials
- Rackable tubes for convenient storage and handling
- Medical grade plastics for assured quality
- Gamma Irradiated for guaranteed sterility



Amies Charcoal

Ideal for prolonged transport of clinically significant samples including fastidious and non-fastidious organisms, with up to 48 hours recovery.

Amies Clear

Recommended for microscope work and general clinical swabbing of fastidious, non-fastidious and anaerobic microorganisms, with up to 24 hours recovery.

Amies Liquid

1ml of liquid Amies medium with flocked swab, offering up to 48 hours recovery for aerobic, anaerobic and fastidious organisms. Ideal for high throughput automated laboratory systems & recommended for general sampling including nasopharyngeal, endocervical and wound sampling.

For further details please visit www.tscswabs.co.uk or call +44 (0)1706 620600

15/04/2014_Issue2

