



## Product Specification Sheet

### **TS/164**

Non-Sterile 2500 ml Empty Sampling Container

Product Code	TS/164
Description	Non-Sterile 2500 ml Empty Sampling Container
Intended Use	Water testing or Urine collection over 24 hours for chemical tests
Components	HDPE Bottle HDPE Cap Label
Contents	Not applicable as supplied empty.
Method of Sterilisation	N/A - supplied unsterile
Shelf Life	N/A
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Precautions and Limitations	The bottle is made of PTFE and is only suitable for water and certain solvents. Please check compatability before use
Pack Size	100 in a box
Other	Quality control testing details on final product:1) Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 5 Date : 06.06.2016

ProductSpec\_TS164