



Product Specification Sheet

TS/18-H500

3 inch polystyrene shaft swab with polyester tip in 3 inch tube.

Product Code	TS/18-H500
Description	3 inch polystyrene shaft swab with polyester tip in 3 inch tube.
Intended Use	Ideal for general swabbing where a small shaft is required.
Components	Polypropylene Tube, HDPE cap, polystyrene/Polyester tip swab, Label.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Using dry swabs only gives limited time of recovery
Pack Size	500 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 0 Date: 07.08.2026

ProductSpec_TS18-H500