



Product Specification Sheet

BAG 1

100 Small Minigrip bags, sealed inside 2 Small Minigrip bags - both outer bags tamper sealed with Irradiation indicators. 10 packs (1000) per box

| | |
|-----------------------------|---|
| Product Code | BAG 1 |
| Description | 100 Small Minigrip bags, sealed inside 2 Small Minigrip bags - both outer bags tamper sealed with Irradiation indicators. 10 packs (1000) per box |
| Intended Use | Sample throughout. General laboratory use only |
| Components | Small Minigrrips (220mm x 130mm), Irradiation sterility indicators x20. Recyclable Food grade LDPE |
| Contents | Not applicable. |
| Method of Sterilisation | Sterilised by gamma irradiation at dosage of 25-53.4kGy |
| Shelf Life | 24 months from date of manufacture |
| Storage Conditions | 4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging |
| Precautions and Limitations | Not applicable. |
| Pack Size | 1000 in a box |
| Other | Quality control testing details on final product: i) Validation of irradiation sterility. |

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 20.10.2020
ProductSpec_BAG1