



Product Specification Sheet

TS/15-5B

Poultry boot swab kit. 5 pairs premoistened in MRD all in resealable outer bag for sample transportation (includes overboots & gloves)

Product Code	TS/15-5B
Description	Poultry boot swab kit. 5 pairs premoistened in MRD all in resealable outer bag for sample transportation (includes overboots & gloves)
Intended Use	Poultry house litter testing for Salmonellae
Components	10 boot swabs – 2 x Plastic overboots - 2 x Minigrip bag - 2 x gloves - 2xMedical grade bag
Contents	De-ionised Water - 1 Litre, Peptone – 1g, Sodium Chloride – 8.5g, pH – 7.0 ± 0.2.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	14 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	1. Peptone Saline Diluent (PSD) has been modified to enhance enrichment of target organisms. Further selective steps may be required. 2. Due to increased material absorbance over time, boots may vary in appearance as the liquid media settles into the deeper fibres leading to coloured precipitate on the visible surface. This is a known characteristic and does not affect performance or sterility. If apparent, gentle massage the boots whilst still sealed to redistribute the media before use.
Pack Size	20 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility ii) Organism Recovery (Salmonella typhimurium).

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 2 Date: 03.09.2024

ProductSpec_TS15-5B