



Product Specification Sheet

TS/15-ZA

Sterile 1000 ml Water Sampling Container with Sodium Thiosulfate (180 mg)

| | |
|-----------------------------|---|
| Product Code | TS/15-ZA |
| Description | Sterile 1000 ml Water Sampling Container with Sodium Thiosulfate (180 mg) |
| Intended Use | Sampling of potable and swimming pool water samples for microbiological analysis. Sodium thiosulphate neutralises chlorine (up to 50mg/L) |
| Components | HDPE Bottle, HDPE Cap, Label. |
| Contents | Sodium thiosulphate- 180mg. |
| Method of Sterilisation | Sterilised by gamma irradiation at dosage of 25-53.4kGy |
| Shelf Life | 12 months from date of manufacture |
| Storage Conditions | 4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging |
| Precautions and Limitations | Not applicable. |
| Pack Size | 100 in a box |
| Other | Quality control testing details on final product: i) Validation of irradiation sterility. |

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 3 Date: 24.04.2020

ProductSpec_TS15-ZA