



Product Specification Sheet

TS/8-A

Sterile Woodshaft/cotton swab, 1 in peel pouch

Product Code	TS/8-A
Description	Sterile Woodshaft/cotton swab, 1 in peel pouch
Intended Use	Culture swab & general laboratory use
Components	Woodshaft/Cotton swab, Peel pouch film.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Certain viruses and Micro-organisms are sensitive to cotton and maybe stressed or not even recover
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 5 Date: 24.04.2020
ProductSpec_TS8-A