



# Technical Service Consultants Ltd

CE and UKCA marked Medical Devices  
Hygiene Monitoring  
Pathogen Testing  
Laboratory Consumables



# Index



## Price List Navigation:

Click on images below to visit the required section.



**Probact Transport Swabs**



**Sterile Dry Swabs**



**Specialised Swabs & Components**



**Diagnostic Swabs**



**Protect Microorganism Preservation System**



**Laboratory Accessories**



**Enviroscreen Range**



**Containers**



## Technical Service Consultants Ltd

CE and UKCA marked Medical Devices  
Hygiene Monitoring  
Pathogen Testing  
Laboratory Consumables



# Proback Transport Swabs

Complete sampling kit for timed transport of all types of Microorganisms to the laboratory. Amies Clear medium gives up to 24 hours survival and is ideal for microscopy work. Amies Charcoal medium gives up to 48 hours for transport since the activated charcoal component absorbs toxic metabolites away from the sample.

Amies Liquid medium is ideal for maintaining the viability of aerobic, anaerobic and fastidious bacteria for up to 48 hours.

A choice of swab buds including synthetic and natural materials allows the user to tailor cost and performance. Gamma irradiated to SAL 10<sup>-6</sup> for guaranteed sterility.

Please check the TSC website for individual product information including guidance for use, applications and limitations.

## AMIES CLEAR - Shelf Life 2 Years


M40-A2 Compliant & Microscope work. Mini size available. Up to 24 hours transport.



Product Code	Swab Tip Type	Swab Shaft	Format	Packaging	Breakpoint	Breaks	Quantity	Price
TS/5-1 <sup>+</sup>	Cotton	Tight Grain Wood	Swab & Media Tube	Peel Pouch	No	Yes	100	£91.01
TS/5-3 <sup>+</sup>	Polyester	Polystyrene	Swab & Media Tube	Peel Pouch	No	Yes	100	£96.65
TS/5-5 <sup>+</sup>	Viscose Mini	Aluminium	Swab & Media Tube	Peel Pouch	No	No	100	£138.07
TS/5-9 <sup>+</sup>	Viscose	Polystyrene	Swab & Media Tube	Peel Pouch	45mm	Yes	100	£95.71
TS/5-17 <sup>+</sup>	Viscose	Polystyrene	Swab & Media Tube	Peel Pouch	No	No	100	£91.96
TS/5-17A <sup>+</sup>	Cotton	3 Inch Polystyrene	Mini Swab & Tube	Peel Pouch	No	No	100	£124.53

## AMIES CHARCOAL - Shelf Life 2 Years


M40-A2 Compliant & extended recovery. Up to 48 hours transport.



Product Code	Swab Tip Type	Swab Shaft	Format	Packaging	Breakpoint	Breaks	Quantity	Price
TS/5-2 <sup>+</sup>	Cotton	Tight Grain Wood	Swab & Media Tube	Peel Pouch	No	Yes	100	£91.01
TS/5-4 <sup>+</sup>	Polyester	Polystyrene	Swab & Media Tube	Peel Pouch	No	Yes	100	£96.65
TS/5-6 <sup>+</sup>	Viscose Mini	Aluminium	Swab & Media Tube	Peel Pouch	No	No	100	£138.07
TS/5-8 <sup>+</sup>	Cotton Mini	Treated Paper	Swab & Media Tube	Peel Pouch	No	Yes	100	£126.74
TS/5-10 <sup>+</sup>	Viscose	Polystyrene	Swab & Media Tube	Peel Pouch	45mm	Yes	100	£95.71
TS/5-18 <sup>+</sup>	Viscose	Polystyrene	Swab & Media Tube	Peel Pouch	No	No	100	£91.96
TS/5-18A <sup>+</sup>	Cotton	3 Inch Polystyrene	Swab & Media Tube	Peel Pouch	No	No	100	£124.53

## CARY BLAIR & STUARTS - Shelf Life 2 Years

For Faecal specimens or *Enterobacteriaceae* transportation.




Product Code	Swab Tip Type	Swab Shaft	Format	Packaging	Breakpoint	Breaks	Quantity	Price
TS/5-21 <sup>+</sup>	Viscose	Polystyrene	Swab & Cary Blair Tube	Peel Pouch	No	No	MOQ Applied*	£1,060.83

\*Minimum order quantity applies - 10 x 100

# Liquid Applications

## LIQUID TRANSPORT - From Shelf Life 12 months

1ml for Automated processing of viruses and bacteria.



Product Code	Swab Tip Type	Swab Shaft	Format	Packaging	Breakpoint	Breaks	Quantity	Price
TS/6-17 <sup>+</sup>	Foam Swab	Polystyrene	Swab & Media Tube	Peel Pouch	79mm	Yes	64	£459.14
TS/5-34A <sup>+</sup>	Polyester	Polystyrene	Swab & Media Tube	Peel Pouch	79mm	Yes	125	£178.63

+ Products marked with this symbol are CE marked to the appropriate directive  
Breaks = Swab bent one or more times to at or beyond 90°

## Sterile Dry Swabs

A wide choice of different diameter swab heads and shaft materials combined with a range of synthetic or natural swab materials allows the user to optimise conditions for both sample recovery and cost. Tamper evident label ensures product security and inclusion of tube provides secure transport after sampling.

Alternative peel pouch configuration gives convenience and cost savings where single use is required. Sterilised by Gamma irradiation to ensure a SAL of at least 10<sup>-6</sup>.

Please check the TSC website for individual product information including guidance for use, applications and limitations.



### POLYESTER SWABS - Shelf Life 2 Years

Synthetic material that provides excellent absorbency and release and suitable for all applications and preferred for molecular techniques.

Product Code	Swab Tip Type	Swab Shaft	Format	Packaging	Breakpoint	Breaks	Quantity	Price
TS/19-G250 <sup>+</sup>	Standard	Polystyrene	Single Swab	Tube	No	Yes	250	£120.41
TS/19-J <sup>+</sup>	Standard	Polystyrene	Single Swab	Peel Pouch	No	Yes	1000	£269.53
TS/19-M250	Standard	Polystyrene	Single Swab - Low ATP	Tube	No	Yes	250	£137.14

### VISCOSE SWABS - Shelf Life 2 Years

Modified cellulose from wood pulp. All round usage with good absorbency and release properties. Suitable for all applications.

Product Code	Swab Tip Type	Swab Shaft	Format	Packaging	Breakpoint	Breaks	Quantity	Price
TS/19-P250 <sup>+</sup>	Standard	Polystyrene	Single Swab	Tube	45mm	Yes	250	£109.35
TS/19-S <sup>+</sup>	Standard	Polystyrene	Single Swab	Peel Pouch	45mm	Yes	1000	£192.49
TS/17-A250 <sup>+</sup>	Standard	Polypropylene	Single Swab	Tube	No	No	250	£98.57
TS/17-B <sup>+</sup>	Standard	Polypropylene	Single Swab	Peel Pouch	No	No	1000	£174.51
TS/19-A250 <sup>+</sup>	Standard	Polystyrene	Single Swab	Tube	No	Yes	250	£106.53
TS/19-B <sup>+</sup>	Standard	Polystyrene	Single Swab	Peel Pouch	No	Yes	1000	£183.91
TS/14-C <sup>+</sup>	Mini	Aluminium	Single Swab	Tube	No	No	100	£126.02

### MINI TIP SWABS - Shelf Life 2 Years

Smaller bud for intricate applications or sampling.

Product Code	Swab Tip Type	Swab Shaft	Format	Packaging	Breakpoint	Breaks	Quantity	Price
TS/14-C <sup>+</sup>	Viscose	Aluminium	Single Swab	Tube	No	No	100	£126.02
TS/16-A250 <sup>+</sup>	Cotton	Treated Paper	Single Swab	Tube	No	Yes	250	£142.08
TS/16-B <sup>+</sup>	Cotton	Treated Paper	Single Swab	Peel Pouch	No	Yes	1000	£305.70
TS/16-D	Cotton	Treated Paper	100 Swabs	Minigrip Bag	No	Yes	1000	£162.25
TS/16-S <sup>+</sup>	Cotton	Treated Paper	Single Swab	Peel Pouch	45mm	Yes	1000	£232.88

+ Products marked with this symbol are CE marked to the appropriate directive  
Breaks = Swab bent one or more times to at or beyond 90°





## Sterile Dry Swabs

### COTTON SWABS - Shelf Life 2 Years

Natural material with good absorbency and acceptable release properties. Recommended for general applications and not suitable for certain specimens.



Product Code	Swab Tip Type	Swab Shaft	Format	Packaging	Breakpoint	Breaks	Quantity	Price
TS/6-A250 <sup>+</sup>	Standard	Tight Grain Wood	Single Swab	Tube	No	Yes	250	£91.10
TS/8-A <sup>+</sup>	Standard	Tight Grain Wood	Single Swab	Peel Pouch	No	Yes	1000	£151.13
TS/9-A	Standard	Tight Grain Wood	Packed in 25's	Peel Pouch	No	Yes	1000	£109.95
TS/9-B <sup>+</sup>	Standard	Tight Grain Wood	Packed in 5's	Peel Pouch	No	Yes	1000	£118.67
TS/9-C <sup>+</sup>	Standard	Tight Grain Wood	Packed in 2's	Peel Pouch	No	Yes	1000	£121.79
TS/9-D	Standard	Tight Grain Wood	Packed in 10's	Peel Pouch	No	Yes	1000	£117.60
TS/16-A250 <sup>+</sup>	Mini	Treated Paper	Single Swab	Tube	No	Yes	250	£142.08
TS/16-B <sup>+</sup>	Mini	Treated Paper	Single Swab	Peel Pouch	No	Yes	1000	£305.70
TS/16-S <sup>+</sup>	Mini	Treated Paper	Single Swab	Peel Pouch	45mm	Yes	1000	£232.88
TS/18-B500 <sup>+</sup>	Standard	Polypropylene	3 inch Swab In Tube	Tube	No	No	500	£221.01

## Specialised Swabs & Components

### BULK SWAB PACKS - Shelf Life 2 Years

Breakable shaft using multiple bends or scissors for custom sample transfer.



Product Code	Swab Tip Type	Swab Shaft	Format	Packaging	Breakpoint	Breaks	Quantity	Price
TS/10	Standard	Tight Grain Wood	100 Swabs in Mini Grip Bag	Mini Grip Bag	No	Yes	5000	£185.80
TS/10-B	Standard	Tight Grain Wood	100 Swabs	Paper Bag	No	Yes	8500	£232.78
TS/11	No Tip	Tight Grain Wood	100 Shafts Only	Minigrip Bag	No	Yes	5000	£183.75
TS/16-D	Mini	Treated Paper	100 Swabs	Minigrip Bag	No	Yes	1000	£162.25
TS/18	Standard	Polypropylene	3 inch Swabs, 100 in a Minigrip Bag	Minigrip Bag	No	Yes	5000	£199.04

+ Products marked with this symbol are CE marked to the appropriate directive  
 Breaks = Swab bent one or more times to at or beyond 90°



## Diagnostic Swabs



### SwabSURE Colour Change Pathogen Detection Kits - Shelf Life 1 Year

Screening kits for a number of common pathogens, comes complete with resealable tube and scored premoist swabs.

SS-L01	L.monocytogenes. For direct detection of pathogenic L.monocytogenes by chromogenic blue colour. Other Listeria/commensals remain colourless. Liquid media in Tube plus Premoist 65mm Breakpoint Neutralising Foam Tip Swab.	Inner Box and PreMoist Swabs	81 Tests	<b>£800.28</b>
SS-L02	Listeria Species. For direct detection of Listeria spp. by black colour. Commensals remain yellow. Liquid media in Tube plus Premoist 45mm Breakpoint Neutralising Viscose Swab.	Inner Box and PreMoist Swabs	100 Tests	<b>£300.79</b>
SS-L03	Salmonella. For direct detection of Salmonella spp. by black colour. Commensals remain purple. Liquid media in Tube plus Premoist 45mm Breakpoint Neutralising Viscose Swab.	Inner Box and PreMoist Swabs	100 Tests	<b>£300.79</b>
SS-L04	Ecoli Coliform. For direct detection of E.coli by chromogenic blue colour. Coliforms turn pink whilst Commensals remain colourless. Liquid media in Tube plus Premoist 45mm Breakpoint Neutralising Viscose Swab.	Inner Box and PreMoist Swabs	100 Tests	<b>£300.79</b>

## Sterile Disposable Inoculating Loops, Needles, Spreaders & Pipettes


Laboratory consumables for accurate and easy throughput of samples. Loops are available in a choice of firm and soft materials and a range of volume sizes to suit user

preference, sample contamination levels or to comply with recommended techniques. Packaging options allow for cost effective high volume use in resealable bags or packaged

individually in a peel pouch for aseptic techniques. Sterilised by Gamma irradiation to ensure a SAL of at least 10-6.

### DISPOSABLE INOCULATING LOOPS - Shelf Life 3 Years

Smooth and quick sample transfer for increased work throughput.



TS/1-A	Inoculating Loops-10µl	White, Soft	20 per bag	1000	<b>£92.19</b>
TS/1-C	Inoculating Loops-10µl	White, Soft	Single in Peel Pouch	1000	<b>£181.24</b>
TS/2-A	Inoculating Loops-10µl	Green, Firm	20 per bag	1000	<b>£92.19</b>
TS/2-C	Inoculating Loops-10µl	Green, Firm	Single in Peel Pouch	1000	<b>£181.24</b>
TS/30-D	Inoculating Loops-5µl	Red, Soft	20 per bag	1000	<b>£92.19</b>
TS/30-G	Inoculating Loops-5µl	Red, Soft	Single in Peel Pouch	1000	<b>£181.24</b>
TS/4-A	Inoculating Loops-1µl	Clear, Firm	20 per bag	1000	<b>£92.19</b>
TS/4-C	Inoculating Loops-1µl	Clear, Firm	Single in Peel Pouch	1000	<b>£181.24</b>
TS/4-D	Inoculating Loops-1µl	Yellow, Soft	20 per bag	1000	<b>£92.19</b>

### STERILE DISPOSABLE INOCULATING NEEDLES - Shelf Life 2 Years


Ideal for extracting Protect Beads or Stab inoculation technique.



TS/3-A	Inoculating Needle	White, Flexible	20 per bag	1000	<b>£92.19</b>
TS/3-C	Inoculating Needle	White, Flexible	Single in Peel Pouch	1000	<b>£181.24</b>

### STERILE DISPOSABLE SPREADERS - Shelf Life 3 Years

Smooth spreading for quick coverage giving increased work throughput.



TS/30-A	"L" Shaped Spreaders	Clear, Firm	5 per bag	1000	<b>£160.35</b>
TS/30-B	"L" Shaped Spreaders	Clear, Firm	Single in Peel Pouch	1000	<b>£231.78</b>
TS/30-C	"L" Shaped Spreaders	Clear, Firm	25 per bag	1000	<b>£125.45</b>
TS/30-ES	"L" Shaped Spreaders	Blue, Soft	25 per bag	1000	<b>£125.45</b>

# Protect Microorganism Preservation System

Cryopreservation system for long term storage of microorganisms including anaerobes, yeast and moulds, dairy and pharmaceutical strains.

Sterilised by Gamma irradiation to ensure a SAL of at least 10<sup>-6</sup>. Customised accessories available.



## PROTECT MULTIPURPOSE - Shelf Life 2 Years

Ideal for Quality Control keeping organism identical throughout. Useful in seed vessel or fermentation applications.

TS/80-MX	Mixed coloured caps & beads	Moisture Resistant Freezer Box		80	<b>£105.78</b>
TS/80-YE	Yellow caps & beads	Moisture Resistant Freezer Box		80	<b>£105.78</b>
TS/80-GN	Green caps & beads	Moisture Resistant Freezer Box		80	<b>£105.78</b>
TS/80-RD	Red caps & beads	Moisture Resistant Freezer Box		80	<b>£105.78</b>
TS/80-WH	White caps & beads	Moisture Resistant Freezer Box		80	<b>£105.78</b>
TS/80-BL	Blue caps & beads	Moisture Resistant Freezer Box		80	<b>£105.78</b>
TS/71-MX	Mixed coloured caps & beads	Polypropylene Tray		100	<b>£164.69</b>
TS/71-YE	Yellow caps & beads	Polypropylene Tray		100	<b>£164.69</b>
TS/71-GN	Green caps & beads	Polypropylene Tray		100	<b>£164.69</b>
TS/71-RD	Red caps & beads	Polypropylene Tray		100	<b>£164.69</b>
TS/71-WH	White caps & beads	Polypropylene Tray		100	<b>£164.69</b>
TS/71-BL	Blue caps & beads	Polypropylene Tray		100	<b>£164.69</b>
TS/3-A	Inoculating Needle	White, Flexible	20 per bag	1000	<b>£92.19</b>
TS/30-D	Inoculating Loops-5µl	Red, Soft	20 per bag	1000	<b>£92.19</b>

**Don't forget your Cryoblock (TS/70-H)**



## PROTECT PLUS - Shelf Life 2 Years

Includes Loops and Needles.

TS/71-MX PLUS	Mixed coloured caps & beads	Polypropylene Tray	Loops & Needles	100	<b>£168.36</b>
TS/71-YE PLUS	Yellow caps & beads	Polypropylene Tray	Loops & Needles	100	<b>£168.36</b>
TS/71-GN PLUS	Green caps & beads	Polypropylene Tray	Loops & Needles	100	<b>£168.36</b>
TS/71-RD PLUS	Red caps & beads	Polypropylene Tray	Loops & Needles	100	<b>£168.36</b>
TS/71-WH PLUS	White caps & beads	Polypropylene Tray	Loops & Needles	100	<b>£168.36</b>
TS/71-BL PLUS	Blue caps & beads	Polypropylene Tray	Loops & Needles	100	<b>£168.36</b>

## PROTECT SELECT - Shelf Life 2 Years

Tailored formulations for delicate or specialised organism applications.



TS/73-AN 25	Protect Select - Anaerobes	Moisture Resistant Freezer Box	Orange cap & beads	25	<b>£44.10</b>
TS/73-DY 25	Protect Select - Dairy	Moisture Resistant Freezer Box	Lilac cap & beads	25	<b>£41.75</b>
TS/73-MF 25	Protect Select - Meat-Free	Moisture Resistant Freezer Box	Clear cap & black beads	25	<b>£41.75</b>
TS/73-YM 25	Protect Select - Yeasts & Moulds	Moisture Resistant Freezer Box	Grey cap & beads	25	<b>£44.10</b>



## PROTECT ACCESSORIES

Consumables for use with Protect.

TS/70-H	Cryoblock	NEW DESIGN: 6 wells in Aluminium/ Rubber surround		Each	<b>£77.49</b>
TS/30-D	Sterile Inoculating Loops 5µl	Red, Soft, 20 per Bag	Minigrip bag	1000	<b>£92.19</b>
TS/30-G	Sterile Inoculating Loops 5µl	Red, Soft, Singles	Peel Pouch	1000	<b>£181.24</b>
TS/3-A	Sterile Inoculation Needles	Clear, Firm, 20 per Bag	Minigrip bag	1000	<b>£92.19</b>
TS/3-C	Sterile Inoculation Needles	Clear, Firm, Single	Peel Pouch	1000	<b>£181.24</b>

## Laboratory Accessories


Consumables and Items for sanitation.



### STERILE BAGS

For aseptic or sterile applications.

Bag 1	125mm x 215mm Sterile	Small Minigrip Bag		1000	<b>£280.65</b>
Bag 3	230mm x 350mm Sterile	Large Minigrip Bag		1000	<b>£296.21</b>





# Enviroscreen Range

A range of environmental monitoring products offered in a variety of different formats to optimise hygiene monitoring protocols.

Sterilised by Gamma irradiation to ensure a SAL of at least 10<sup>-6</sup>.

A choice of neutralising products are available for use in areas where traces of disinfectants may hide contamination.



## SWAB SAMPLING KIT - Shelf Life 2 Years

All in one crevice and surface sampling with disinfectant neutralisers included.



TS/5-40	5ml Neutralising Buffer	Viscose Tip, Blue Polystyrene Breakpoint Shaft	Peel Pouch	100	£98.84
TS/5-42	10ml Neutralising Buffer	Polyester Tip, Polystyrene Shaft	Peel Pouch	100	£132.74

**Don't forget - when sampling larger areas, use a Sponge Kit**

## PREMOIST SWAB KIT - Shelf Life 2 Years

Pre wetted swab for convenience in resealable tube for transport.



TS/6-61:250	Premoistened in TSC Neutralising Buffer	Viscose Tip, Polystyrene Shaft	Tube	250	£112.63
TS/6-61BP250	Premoistened in TSC Neutralising Buffer	Viscose Tip, Blue Polystyrene Breakpoint Shaft	Tube	250	£114.20
TS/6-62:100	Premoistened in Deionised Water	Viscose Tip, Polystyrene Shaft	Tube	100	£41.16
TS/6-63BP100	Premoistened in TSC Neutralising Buffer	Foam Tip, Polystyrene Breakpoint Shaft	Tube	100	£71.07

## MAXIMUM RECOVERY DILUENT KIT - Shelf Life 2 Years

Swab with no neutralisers included.



TS/5-31	Liquid MRD Medium	Viscose Tip, Polystyrene Shaft	Tube in Peel Pouch	100	£117.56
---------	-------------------	--------------------------------	--------------------	-----	---------

## ENVIROSCREEN DRY SWABS - Shelf Life 2 Years

Non-inhibitory swab supplied in resealable tube for sample protection.



TS/6-H250	Wood Shaft, Cotton Tip	Hygiene Woodstick Cotton Tip Swabs	Tube	250	£91.10
TS/15-CP	12 inch Plastic Shaft	Blue Jumbo Foam Swab Tip	Tube	100	£207.34
TS/17-H250	Polypropylene Shaft, Viscose Tip	Polypropylene Shaft Viscose Tip	Tube	250	£98.57
TS/19-H250	Blue Polystyrene Breakpoint Shaft, Viscose Tip	Blue Polystyrene Breakpoint Shaft, Viscose Tip	Tube	250	£109.96
TS/19-T	Blue Polystyrene Breakpoint Shaft, Viscose Tip	Blue Polystyrene Breakpoint Shaft, Viscose Tip	Peel Pouch	1000	£252.54
TS/19-M250	Ultra Clean Swab for Bioluminescence	Polystyrene Shaft, Polyester Tip	Tube	250	£137.14

## SAMPLING TEMPLATES - Shelf Life 2 Years

Ensures sampling technique to recommended & approved standards.



TS/15-T40	Sterile Plastic Template	10 x 10cm <sup>2</sup> Blue Plastic	Individually Wrapped	40	£70.78
-----------	--------------------------	-------------------------------------	----------------------	----	--------



## ENVIROSCREEN SPONGE POUCHES - Shelf Life 14 Months

All-in-one method with no sponge transfer minimising potential sample contamination.

TS/15-B:NaCl	50cm <sup>2</sup> Blue Sponge	Dosed in Saline (0.9% NaCl)	Easy Open Stomacher Pouch	200	£336.67
TS/15-SA	EnviroStik - Blue Sponge	Dosed in 0.9% NaCl, Minigrip Bag	Easy Open Stomacher Pouch	125	£239.61
TS/15-B:PSD	50cm <sup>2</sup> Blue Sponge	Dosed in Peptone Saline Diluent (MRD)	Easy Open Stomacher Pouch	200	£313.21
TS/15-BPW	50cm <sup>2</sup> Blue Sponge	Dosed in Buffered Peptone Water	Easy Open Stomacher Pouch	200	£327.89
TS/15-B	50cm <sup>2</sup> Blue Sponge	Dosed in Neutralising Buffer	Easy Open Stomacher Pouch	200	£342.63
TS/15-SH	EnviroStik - Blue Sponge	Dosed in Neutralising Buffer, Minigrip Bag	Easy Open Stomacher Pouch	125	£254.11
TS/15-YCS	100cm <sup>2</sup> Yellow Polyurethane Sponge	Dosed in Peptone Saline Diluent	Easy Open Stomacher Pouch	75	£126.25



## ENVIROSCREEN SPONGE CONTAINERS - Shelf Life 1 Year

Alternative All in one method with no sponge transfer in stackable format.

TS/15-A-NaCl95	50cm <sup>2</sup> Blue Sponge	Dosed in Saline (0.9% NaCl)	100ml Plastic Container	MOQ Applied*	£372.07
TS/15-A-95	50cm <sup>2</sup> Blue Sponge	Dosed in Neutralising Buffer	100mm Plastic Container	MOQ Applied*	£402.23

\*Minimum order quantity applies - 5 boxes of 95



## POULTRY SAMPLING BOOT SWABS - Shelf Life 14 months

Environmental sampling & collection kits for the presence of Salmonellae & Campylobacter spp.

TS/15-2A	Poultry Boot Swab 2 Pairs	Dosed in MRD, Gloves + Bag	2x Easy Open Stomacher Bags	30	£65.69
TS/15-2B	Poultry Boot Swab 2 Pairs	Dosed in MRD, Gloves, Bag, Overboots	2x Easy Open Stomacher Bags	25	£104.63
TS/15-5A	Poultry Boot Swab 5 Pairs	Dosed in MRD, Gloves + Bag	2x Easy Open Stomacher Bags	20	£84.06
TS/15-5B	Poultry Boot Swab 5 Pairs	Dosed in MRD, Gloves, Bag, Overboots	2x Easy Open Stomacher Bags	20	£91.47



## MEAT CARCASS SAMPLING KITS - Shelf Life 14 months

Direct food Contact sampling for pathogenic organisms including Salmonellae.

TS/15-B:NaCl	50cm <sup>2</sup> Blue Sponge	Dosed in Saline (0.9% NaCl)	Easy Open Stomacher Pouch	200	£336.67
TS/15-SA	Envirostik Blue Sponge on Stick	Dosed in 0.9% NaCl	Easy Open Stomacher Pouch	125	£239.61
TS/15-B:PSD	50cm <sup>2</sup> Blue Sponge	Dosed in Peptone Saline Diluent (MRD)	Easy Open Stomacher Pouch	200	£313.21
TS/15-BPW	50cm <sup>2</sup> Blue Sponge	Dosed in Buffered Peptone Water	Easy Open Stomacher Pouch	200	£327.89
TS/15-BPS	EnviroStik - Blue Sponge	Dosed in Buffered Peptone Water	Easy Open Stomacher Pouch	125	£249.60

# Containers

A range of food safe polyethylene sterile and non sterile containers for sampling from a variety of environments. Available with the option of sodium thiosulphate included for chlorine removal from potable water samples.



**LABELLED CONTAINERS - Unsterile for Urine Collection, easy screw and leakproof tested**  
 Suitable for Urine collection, unsterile and Graduated for information. Now with updated label.

TS/164 <sup>+</sup>	2500ml	2.5 Litre 24 Hour Labelled Urine Container	Polyethylene	100	£177.81
---------------------	--------	--	--------------	-----	---------



**WATER SAMPLING CONTAINERS - From Shelf Life 1 Year**  
 Sterile, with & without tamper evident, final volume graduated. Available with Neutraliser for removal of free & combined residual chlorine (as per SCA recommendations).

TS/15-F40	500ml Tamper Evident	Empty (2 year expiry)	n/a	40	£75.87
TS/15-D40	500ml Tamper Evident	60mg Sodium Thiosulphate	Neutralise 28mg Chlorine/Litre	40	£87.46
TS/15-DA40	500ml Tamper Evident	90mg Sodium Thiosulphate	Neutralise 50mg Chlorine/Litre	40	£117.33
TS/15-DL40	500ml Tamper Evident	9mg Sodium Thiosulphate	Neutralise 5mg Chlorine/Litre	40	£114.92
TS/15-Z	1000ml	120mg Sodium Thiosulphate	Neutralises 33mg Chlorine/Litre	100	£224.06
TS/15-ZA	1000ml	180mg Sodium Thiosulphate	Neutralise 50mg Chlorine/Litre	100	£233.64

+ Products marked with this symbol are CE marked to the appropriate directive

# Plastic Injection Moulding

We are pleased to offer quotations for any Injection Moulding requirements you may have. Goodbrand Plastics are manufacturers of Plastic Injection Mouldings, catering for large to small batches. Components are produced for a varied range of industries including Medical, Construction, Highways, Security, Automotive, Leisure & Electronics. In house assembly is available if required.

Please contact Goodbrand Plastics on **01270 753006** or e-mail **graham@goodbrandplastics.co.uk**

For more information please visit our website at: **www.goodbrandplastics.co.uk**

### Interest Charges

A finance charge of 1.5% per month is charged on all balances that remain unpaid beyond payment terms.

### Summary of terms and conditions

Full terms and conditions can be downloaded from our website

- 1 Prices are quoted Pound Sterling and are net, ex-works and subject to alteration without notice.
- 2 VAT will be charged at the rate prevailing at date of invoice.
- 3 Invoices due for payment 30 days from date of issue.
- 4 Discounts are available for bulk or standing orders.
- 5a Delivery to mainland UK on a next day service at £40.00 up to 20kg. A sliding scale of costs operates on higher weights. Requests for bulk or timed deliveries will be treated individually and costs agreed prior to dispatch.
- 5b Northern Ireland carriage please enquire. A sliding scale of costs operates on higher weights.
- 5c Supplement will be incurred for Western Isles.
- 6 All goods will remain the property of TSC Ltd. until the price has been paid in full. All materials are sold for laboratory diagnostic or allied use only and are not intended for administration as a drug.



## Technical Service Consultants Ltd

CE and UKCA marked Medical Devices  
Hygiene Monitoring  
Pathogen Testing  
Laboratory Consumables