



Product Specification Sheet

TS/19-H

Swab with a blue breakable polystyrene shaft and pure viscose tip in a polypropylene tube

Product Code	TS/19-H
Description	Swab with a blue breakable polystyrene shaft and pure viscose tip in a polypropylene tube
Intended Use	Culture swab; general Laboratory use - food use
Components	Polypropylene Tube Cap/pure viscose tip/polystyrene shaft swab Non-toxic glue Tube label Box label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42 KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	500 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 18.02.2008

ProductSpec_TS19-H