



## Product Specification Sheet

### **TS/19-T**

Viscose tipped swab with blue break point polystyrene shaft (@45mm) in a peel pouch

Product Code	TS/19-T
Description	Viscose tipped swab with blue break point polystyrene shaft (@45mm) in a peel pouch
Intended Use	Culture swab; general Laboratory use - food use
Components	Viscose tip/polystyrene shaft swab Non-toxic glue Peel pouch material Box label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42 KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 11.02.2010

ProductSpec\_TS19-T