



Product Specification Sheet

TS/16-P

Sterile Cotton mini tip and plasticised paper shaft swab, breakpoint at 45mm. Single in tube.

Product Code	TS/16-P
Description	Sterile Cotton mini tip and plasticised paper shaft swab, breakpoint at 45mm. Single in tube.
Intended Use	Culture swab; recommended for urethral sampling
Components	Swab with cotton mini tip and plasticised paper shaft, breakpoint at 45mm; polypropylene tube; polyethylene cap; Non-toxic glue; Label & film;
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	500 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 10.01.2008

ProductSpec_TS16-P