



Product Specification Sheet

TS/16-S

Sterile dry swab, single in peel pouch. Cottonmini tip and paper shaft swab, breakpoint at 45mm.

Product Code	TS/16-S
Description	Sterile dry swab, single in peel pouch. Cottonmini tip and paper shaft swab, breakpoint at 45mm.
Intended Use	Culture swab; recommended for urethral sampling
Components	Swab with cotton mini tip (tip diameter 1.5 mm) and plasticised paper shaft, breakpoint at 45mm Plastic peel pouch Non-toxic glue Label &film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect fromlight, heat and moisture in tightly closedoriginal container
Pack Size	1000 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 27.12.2012

ProductSpec_TS16-S