



## Product Specification Sheet

### **TS/19-S**

Viscose tip, polystyrene shaft, breakpoint at 45mm, in vacuum packed peel pouch.

Product Code	TS/19-S
Description	Viscose tip, polystyrene shaft, breakpoint at 45mm, in vacuum packed peel pouch.
Intended Use	Culture swab; general Laboratory use
Components	White cap Polypropylene Tube White cap/viscose tip/polystyrene shaft swab (BP 45mm) (tip diameter – 4.5-5.0 mm) Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42K Gy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 27.12.2012

ProductSpec\_TS19-S