



Product Specification Sheet

TS/25

Female & Male Swab pack

Product Code	TS/25
Description	Female & Male Swab pack
Intended Use	Collection and transport of specimens - dual swab pack
Components	Cotton mini tip/plasticised paper shaft(Male); Viscose/Polystyrene (Female); Peel Pouch; Label &film;
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect fromlight, heat and moisture in tightly closedoriginal container
Pack Size	1000 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 3 Date : 27.12.2012

ProductSpec_TS25