



Product Specification Sheet

TS/14-H

Sterile dry swab with reduced size cotton tip and Supreme Twisted wire shaft, single in polypropylene tube

Product Code	TS/14-H
Description	Sterile dry swab with reduced size cotton tip and Supreme Twisted wire shaft, single in polypropylene tube
Intended Use	Collection and transport of specimens
Components	Swab with reduced size cotton tip and twisted wire shaft; mounted in a plastic orange cap; Polypropylene tube; Non-toxic glue; Label & film;
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 23.03.2007

ProductSpec_TS14-H