



Product Specification Sheet

TS/16-A

Sterile dry swab with cotton mini tip and papershaft single in tube.

Product Code	TS/16-A
Description	Sterile dry swab with cotton mini tip and papershaft single in tube.
Intended Use	Culture swab; recommended for urethral sampling
Components	Swab with cotton mini tip (tip diameter 1.5 mm) and plasticised paper shaft, mounted in a polyethylene cap Polypropylene Tube Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	500 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 11.11.2007

ProductSpec_TS16-A