



## Product Specification Sheet

### **TS/16-C**

Sterile dry swab with cotton (charcoal coated) mini tip and paper shaft single in peel pouch.

Product Code	TS/16-C
Description	Sterile dry swab with cotton (charcoal coated) mini tip and paper shaft single in peel pouch.
Intended Use	Culture swab; recommended for urethral sampling
Components	Swab with cotton mini tip and plasticised paper shaft, mounted in a polyethylene cap, polypropylene tube; Activated Charcoal, Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	500 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 11.05.2008

ProductSpec\_TS16-C