



Product Specification Sheet

TS/16-D

Sterile dry swab with cotton tip and plasticised paper shaft, 100 swabs in Minigrip bag.

Product Code	TS/16-D
Description	Sterile dry swab with cotton tip and plasticised paper shaft, 100 swabs in Minigrip bag.
Intended Use	Culture swab; recommended for urethral sampling
Components	Cotton tip/plasticised Paper shaft swab; Minigrip bag; Label & film;
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42K Gy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 27.12.2012

ProductSpec_TS16-D