



Product Specification Sheet

TS/19-J

Swab with a breakable polystyrene shaft & polyester tip in a peel pouch

Product Code	TS/19-J
Description	Swab with a breakable polystyrene shaft & polyester tip in a peel pouch
Intended Use	Culture swab; Patient sampling & general Laboratory use
Components	Polyester tip/polystyrene shaft swab Non-toxic glue Peel pouch material Box label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42 KGy
Shelf Life	24 months from date of manufacture Storage Conditions - 4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 3 Date : 27.12.2012

ProductSpec_TS19-J