



Product Specification Sheet

TS/17-B

Swab with a flexible polypropylene shaft tipped with pure viscose, single in peel pouch

Product Code	TS/17-B
Description	Swab with a flexible polypropylene shaft tipped with pure viscose, single in peel pouch
Intended Use	Culture swab; Patient sampling & general Laboratory use
Components	Polystyrene cap viscose tip/polypropylene shaft swab (tip diameter – 5.15 mm) Peel pouch Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 27.12.2012

ProductSpec_TS17-B