



Product Specification Sheet

TS/17-D

Swab with a flexible polypropylene shaft/cotton tip in a polypropylene tube.

Product Code	TS/17-D
Description	Swab with a flexible polypropylene shaft/cotton tip in a polypropylene tube.
Intended Use	Culture swab; Patient sampling & general Laboratory use
Components	Polypropylene Tube Polyethylene cap cotton tip/polypropylene shaft swab Non-toxic glue Bag label Box label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	500 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 09.05.2007

ProductSpec_TS17-D