



## Product Specification Sheet

### **TS/17-DC**

Swab with a flexible polypropylene shaft and cotton tip in a polypropylene tube. 50 tubes in the sealed outer bag

Product Code	TS/17-DC
Description	Swab with a flexible polypropylene shaft and cotton tip in a polypropylene tube. 50 tubes in the sealed outer bag
Intended Use	Clean Room monitoring. Customised product
Components	White cap Polypropylene Tube cotton tip/polypropylene shaft swab Non-toxic glue Plastic bag Bag label Box label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42 KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	500 (10 packs of 50) in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 09.05.2007

ProductSpec\_TS17-DC