



## Product Specification Sheet

### **TS/18**

3 inch cotton tipped swab with a flexible polypropylene shaft

Product Code	TS/18
Description	3 inch cotton tipped swab with a flexible polypropylene shaft
Intended Use	Ideal for general swabbing where a small shaft is required.
Components	Cotton tip/polypropylene shaft swab (tip diameter – 5.15 mm) Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	5000 in a box (50 packs of 100)
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 18.09.2012

ProductSpec\_TS18