



## Product Specification Sheet

### **TS/10-B**

Sterile swab with a wooden shaft and cottontip, 100 in sealed paper bag

Product Code	TS/10-B
Description	Sterile swab with a wooden shaft and cottontip, 100 in sealed paper bag
Intended Use	Sample throughout. General laboratory use only
Components	Cotton tip/ wooden shaft swab; paper bag; Non-toxic glue; Label & film;
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42K Gy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	8500 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 09.04.2008

ProductSpec\_TS10-B