



Product Specification Sheet

TS/10-NS

Non-Sterile swab with a wooden shaft and cotton tip, 100 in the mini grip bag

Product Code	TS/10-NS
Description	Non-Sterile swab with a wooden shaft and cotton tip, 100 in the mini grip bag
Intended Use	General laboratory use only
Components	Cotton tip/ wooden shaft swab Mini grip bag Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Non-sterile
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	5000 in a box
Other	Quality control testing details on final product:i)

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 14.09.2012

ProductSpec_TS10-NS