



## Product Specification Sheet

### **TS/6-A**

Sterile swab with a wooden shaft and cottontip, in polypropylene outer tube

Product Code	TS/6-A
Description	Sterile swab with a wooden shaft and cottontip, in polypropylene outer tube
Intended Use	Culture swab; general laboratory use
Components	Polypropylene TubeSwab with Blue cap Cotton tip/wooden shaft Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect fromlight, heat and moisture in tightly closedoriginal container
Pack Size	500 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 30.05.2006

ProductSpec\_TS6-A