



## Product Specification Sheet

### **TS/7**

Sterile swab with a wooden shaft and calcium alginate tip, in polypropylene outer tube

Product Code	TS/7
Description	Sterile swab with a wooden shaft and calcium alginate tip, in polypropylene outer tube
Intended Use	General sampling. When used with Calgon Ringers swab head dissolves completely to allow complete sample recovery
Components	White cap Polypropylene Tube Cap calcium alginate tip/woodstick shaft swab Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 16.07.2007

ProductSpec\_TS7