



## Product Specification Sheet

### **TS/8-B**

Sterile swab with a wooden shaft and cotton (charcoal coated) tip, single in peel pouch.

Product Code	TS/8-B
Description	Sterile swab with a wooden shaft and cotton (charcoal coated) tip, single in peel pouch.
Intended Use	Culture swab & general laboratory use
Components	Cotton tip/ wooden shaft swab Activated Charcoal Peel pouch Non-toxic glue Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	1000 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 2 Date : 27.12.2012

ProductSpec\_TS8-B