



Product Specification -, **TS/15-T**, Sterile plastic template 10x10 cm²

Product Code	TS/15-T
Description	Sterile plastic template 10 x 10cm ²
Intended Use	Controlled surface sampling for environmental monitoring (compliant to ISO 18593)
Components	Blue polypropylene, Medical grade Bag, Box Label.
Contents	N/A
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	50 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date :22.09.2010
Spec_TS/15-T

