



## Product Specification Sheet

### **TS/5-40BPW**

5ml of Buffered Peptone Water in tube with bluebreakpoint swab in a peel pouch

Product Code	TS/5-40BPW
Description	5ml of Buffered Peptone Water in tube with bluebreakpoint swab in a peel pouch
Intended Use	Swab for pathogen monitoring in poultry houses Can be used for sampling of non-disinfected surfaces or hand swabbing.
Components	Blue tulip cap Blue viscose/polystyrene breakpoint swab Non-toxic glue Tube Label & film Buffered Peptone Water
Contents	De-ionised Water, Peptone Water Powder, Phosphate Buffers
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism Recovery

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 12.04.2007

ProductSpec\_TS5-40BPW