



## Product Specification Sheet

### **TS/5-62**

Swab for environmental sampling, premoistened in deionised water in tube in a peel pouch

Product Code	TS/5-62
Description	Swab for environmental sampling, premoistened in deionised water in tube in a peel pouch
Intended Use	Environmental sampling. Not suitable for transport. For cleaning validation and hygiene monitoring the samples should be processed on growth media with neutralisers
Components	Blue cap/white viscose tip/polystyrene shaft swab 150 mm Polypropylene Tube, Tube Label Non-toxic glue Base web film and plain film De-ionised water (200 µl)
Contents	De-ionised Water - 1 Litre
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42 KGy
Shelf Life	12 months from the date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility. ii) Product Appearance

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 3 Date : 24.10.2008

ProductSpec\_TS5-62