



Product Specification Sheet

TS/6-H

Sterile swab with a wooden shaft and cottontip, in polypropylene outer tube.
(Hygiene)

Product Code	TS/6-H
Description	Sterile swab with a wooden shaft and cottontip, in polypropylene outer tube. (Hygiene)
Intended Use	Culture swab; general laboratory use
Components	Polypropylene TubeSwab with Blue cap Cotton tip/wooden shaft Non-toxic glue Hygiene Label & film
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	500 in a box
Other	Quality control testing details on final product:i) Validation of irradiation sterility.

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 11.03.2008

ProductSpec_TS6-H