



Product Specification **BAG1**, 100 Small Minigrip bags, sealed inside 2 Small Minigrip bags - both outer bags tamper sealed with Irradiation indicators. 10 packs (1000) per box

Product Code	BAG1
Description	100 Small Minigrip bags, sealed inside 2 Small Minigrip bags - both outer bags tamper sealed with Irradiation indicators. 10 packs (1000) per box
Intended Use	Sample throughout. General laboratory use only
Components	Small Minigrips x 110 (125mm x 125mm); Irradiation Sterility indicators x20;
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from date of manufacture
Storage Conditions	4-25 °C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	1000 in a box
Other	Quality control testing details on final product: i) Validation of irradiation sterility.

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 09.04.2008
ProductSpec_BAG1

