



Product Specification Sheet

TS/126-H

Blown Plastic pipette, Ultra fine tip (50 Dropper). Single in peel pouch

Product Code	TS/126-H
Description	Blown Plastic pipette, Ultra fine tip (50 Dropper). Single in peel pouch
Intended Use	Sample throughout or Protect. General laboratory use only
Components	Clear Plastic (soft) peel pouch
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42K Gy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	500 in a box
Other	Quality control testing details on final product: Validation of irradiation sterility

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 11.03.2008

ProductSpec_TS126-H