



Product Specification **TS/69-MX 10**, Protect Starter Kit - Preservation system for General microorganisms, 5 colours (Protect General use)

Product Code	TS/69-MX 10
Description	Protect Starter Kit - Preservation system for General microorganisms, 5 colours (Protect General use)
Intended Use	Long-term storage of anaerobic microorganisms at low temperatures
Components	Cryovials with same coloured cap and Ceramic beads (20-25); Minigrip Bag; Cryopreservative fluid
Contents	Cryopreservative fluid: Beef extract Peptone Sodium chloride Glycerol (20%) De-ionised Water
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-42KGy
Shelf Life	2 years from the date of manufacture
Storage Conditions	4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	10 in a Bag
Other	Quality control testing details on final product:i) Validation of irradiation sterility.ii) Organism recovery Protect Select - Dairy

This specification is intended for guidance only. You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 1 Date : 23.09.2008

ProductSpec_TS69-MX 10

