



Product Specification Sheet

TS/71 (Covers all products)

Cryovials containing ceramic beads plus a cryopreservative fluid.

Product Code	TS/71 (Covers all products)
Description	Cryovials containing ceramic beads plus a cryopreservative fluid.
Intended Use	Long-term storage of microorganisms at low temperatures: -20°C; -70°C; -80°C; liquid nitrogen.
Components	Cryovials with coloured cap 20-25 Ceramic beads Sterile Loops (Protect PLUS only) Sterile Needles (Protect PLUS only) Polypropylene tray (TS/71 products) Cryopreservative fluid
Contents	Cryopreservative fluid: Beef extract Peptone Sodium chloride Glycerol (20%) De-ionised Water
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 12.5-42K Gy
Shelf Life	24 months from the date of manufacture
Storage Conditions	For unopened vials: 4-25 °C; in dry and ventilated area; protect from light, heat and moisture in tightly closed original container
Pack Size	100 vials in a box
Other	Quality control testing details on final product: Validation of irradiation sterility. Recovery of <i>N. gonorrhoeae</i> Recovery of <i>H. influenzae</i>

This specification is intended for guidance only.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue # : 6 Date : 07.09.2012

ProductSpec_TS/71